

Recommendations for Safe and Long-Lasting Use of Your Device

Dear Customers,

Thank you very much for purchasing our high-quality device. To ensure long-term satisfaction and reliable performance, we would like to draw your attention to several important recommendations:

1. Limiting the Charging Time

Please make sure that the charging time does not exceed 2.5 hours. Longer charging may negatively affect the battery, reduce its lifespan, and in the worst case lead to overheating.

2. Safety During Charging

- Never charge the device under pillows, blankets, or in poorly ventilated areas. Insufficient air circulation during charging may cause overheating and create a fire hazard. Always charge the device on a firm, flat surface with good airflow.

- Do not use powerful or fast-charging adapters (chargers with high output current).

Reason: The battery is designed for a specific charging current. A charger that supplies too much current may cause overheating, accelerate battery wear, or damage the battery.

Use only suitable power adapters with the following parameters: USB power adapter, 5V / 1A.

3. Storage and Transport

For safe storage and transportation, always use an appropriate protective case. It protects the device from scratches, dust, and mechanical damage, ensuring long-term durability.

4. Additional Recommendations

- Avoid drops, strong vibrations, and contact with liquids.
- Do not store the device in humid environments or under direct sunlight.
- Do not leave the device connected to the charger overnight.

By following these simple but important guidelines, you ensure safe operation and significantly extend the lifespan of your device.

Kind regards,

The Biomedis Team